

# Aria 104

## 10" Portable intensive care ventilator

Oxygen driven ventilator with built-in turbine  
for adults, children and newborn

- Touch Screen -



Code: 980104/A

Rev.1- 27/02/2020

### GENERAL DATA

Aria 104 electronic lung ventilator is equipped with turbine and with a TFT 10,4" colour monitor touch screen displaying the curves of pressure, flow, volume, the loops of breathing parameters, the trends and the ventilation parameters. Aria 104 lung ventilator is suitable for ventilation of adult, paediatric and neonatal patients. Aria 104 lung ventilator is equipped with a flow generation system by turbine with separate cooling system granting higher quality and safety standards in patient ventilation.

Aria 104 is equipped with a flow and pressure trigger, it provides the most advanced volume-controlled ventilation modalities VC/VAC, VC/VAC-BABY, pressure-controlled ventilation modalities APCV (BILEVEL ST), APCV-TV, SIMV by Volume and by Pressure, Pressure supported modalities PSV (BILEVEL S), PSV-TV, CPAP, APRV, SIGH, Non-Invasive Ventilation (NIV APCV - NIV PSV), Drug Nebulizer and Manual Ventilation (MAN).

Aria 104 is supplied with back up long lasting batteries and its software can be updated for new modes and last generation ventilation strategies.

### NORMS



The lung ventilator is conform to the essential requirements and it is realized according to the references of the Annex II of 93/42/EEC Medical Devices Directive.

Class and type according to IEC 601-1 Class I Type B

Class according to 93/42 EEC Directive Class IIb

Electromagnetic compatibility (EMC) EN 60601-1-2: 2015 and following

Norms DIR. 93/42/CEE (2007); EN 60601-1 :2006/A1 :2011/A1 :2013; EN 60601-1-2 :2015; IEC 601-1-6:2013; IEC 601-1-8:2012; EN 62304:2006/AC:2008; ISO 10993-1:2009; IEC 62353:2014; EN 60601-2-12:2007; ISO 80601-2-12:2011; ISO 15223-1:2016; DIR. 2011/65/CE; D.Lgs 49/2014; ISO 14971:2012; EN ISO 4135:2001

## ENVIRONMENTAL CONDITIONS

**Operating**                      Relative humidity: 30 - 95% non-condensing  
    Temperature: from -10 to +40°C

**Storage**                        Relative humidity: < 95%  
    Temperature: from -25 to +70°C

## TECHNICAL DATA

**Dimensions (W x H x D)**    290 x 245 x 215 mm

**Weight**                         5,5 Kg

**Electric power supply**       100 - 240Vac / 50 - 60Hz

*Power*                            Max 60 VA

*External power supply*    12 Vdc / 7 A  
*(low voltage)*

*Internal battery*            Battery NiMh 12Vdc - 4.2 Ah

*Internal battery operation*   Max 4 hours

*Battery re-charging time*   About 10 hours

**External electric connections**

- O2 sensor connection
- Flow sensor connection
- CO2 module connection (RS232)
- CPU programming connector (USB 1)
- Data transfer connection: patient data, events, trends (USB 2)
- External alarm/nurse call

**Patient connections**        Male conic connectors 22 mm / Female of 15 mm (according to EN ISO 5356-1:2015 norm)

**Supply pressure (O<sub>2</sub>)**

- Low pressure (max 15 l/min)
- High pressure 280 kPa - 600 kPa / 2.8 - 6 bar / 40 - 86 psi

*Max flow requested (O<sub>2</sub>)*   80 l/min

**IP degree of protection**    IP21

## LUNG VENTILATOR FUNCTIONAL FEATURES

Intended use	Aria 104 is a lung ventilator for use in emergency rooms, transport, intensive care units and with patients affected by respiratory diseases and it is suitable for ventilation of adult, paediatric and neonatal patients.
Operation principle	<ul style="list-style-type: none"> <li>• Time cycled at constant volume</li> <li>• Pressure cycled</li> <li>• Microprocessor controlled flow</li> <li>• Spontaneous breath with integrated valve</li> </ul>
Pressure automatic compensation (altitude)	Automatic compensation of atmospheric pressure on measured pressure: present (max. 5000 mt)
Automatic leaks compensation	Max 60 l/min (NIV APCV, NIV PSV)
Leak % visualization	Present
Visualization of the oxygen consumption calculation	Present
Altitude compensation for oxygen sensor	Present
Respiratory parameters default setting	Present (Neonatal, Paediatric, Adult)
Ventilation modalities	<ul style="list-style-type: none"> <li>• APCV (BILEVEL ST), APCV-TV, PSV (BILEVEL S), PSV-TV (Auto Weaning), VC/VAC, VC/VAC BABY, V-SIMV+PS, P-SIMV+PS, CPAP, APRV</li> <li>• SIGH, NEB (Nebulizer), Apnea BACK-UP (PSV, PSV-TV, CPAP), MAN (Manual).</li> </ul>
Breathing rate VC/VAC	From 4 to 150 bpm
Inspiratory Time / Expiratory Time (maximum, minimum)	<ul style="list-style-type: none"> <li>• Ti min = 0.036sec (minimum inspiratory time)</li> <li>• Ti max = 9.6sec (maximum inspiratory time)</li> <li>• Te min = 0.08sec (minimum expiratory time)</li> <li>• Te max = 10.9sec (maximum expiratory time)</li> </ul>
Breathing rate V-SIMV e P-SIMV	From 1 to 60 bpm
SIMV Inspiratory time	From 0.2 to 5.0 sec.
Tidal volume	<ul style="list-style-type: none"> <li>• From 100 to 3000 ml (Adult)</li> <li>• From 50 to 400 ml (Paediatric)</li> <li>• From 2 to 100 ml (Neonatal)</li> </ul>
I:E ratio	From 1:10 to 4:1

Inspiratory pause	From 0 to 60 % of the inspiratory time
Inspiratory pressure limit	P <sub>insp</sub> : from 2 to 80 cmH <sub>2</sub> O (in function of low and high pressure alarm set)
Inspiratory ramp Slope	1, 2, 3, 4 (acceleration slope) - (4 max. acceleration) (in operative modes by pressure only)
PEEP	From OFF, 2 to 50 cmH <sub>2</sub> O Microprocessor controlled
O <sub>2</sub> concentration	Adjustable from 21 to 100% with electronic integrated mixer.
Trigger detective method	Through sensor (Pressure or Flow)
<i>Pressure trigger ( I )</i>	Pressure adjustable from OFF; -1 to -20 cmH <sub>2</sub> O under PEEP level (step of 1 cmH <sub>2</sub> O)
<i>Flow trigger ( I )</i>	Flow adjustable from OFF; 0.3 to 15 L/min <ul style="list-style-type: none"> <li>• from 0.3 to 1 L/min (step of 0.1 L/min)</li> <li>• from 1 L/min to 2 L/min (step of 0.5 L/min)</li> <li>• from 2 L/min to 15 L/min (step of 1 L/min)</li> </ul>
<i>Trigger E</i>	From 5 to 90 % of the inspiratory peak flow
Inspiratory flow (FLOW)	190 l/min
Flow-by	Automatic
PS (pressure support)	From 2 to 80 cmH <sub>2</sub> O (PSV, V-SIMV+PS, P-SIMV+PS)
SIGH in VC/VAC modality	<ul style="list-style-type: none"> <li>• Interval: 40 ÷ 500 bpm (step 1 bpm)</li> <li>• Amplitude: OFF, 10 ÷ 100% of set Tidal Volume (step 10%)</li> </ul>
CPAP	Pressure: from 3 to 50 cmH <sub>2</sub> O
APRV	<ul style="list-style-type: none"> <li>• Time High and Time Low: from 1 to 200 sec.</li> <li>• Pressure High and Pressure Low: from 3 to 50 cmH<sub>2</sub>O.</li> </ul>

---

Functions	<ul style="list-style-type: none"> <li>• MENU function (SETUP – PATIENT DATA)</li> <li>• Alarms Limits</li> <li>• Graphics visualization (Auto-Range)</li> <li>• INSP Hold - EXP Hold (max 20 sec.)</li> <li>• O<sub>2</sub> 100% control (O<sub>2</sub> to 100% max. 5 min.)</li> <li>• NEB (6 l/min)</li> <li>• MAN (manual ventilation)</li> </ul>
-----------	---

---

NEB	Drug nebulizer: selectable to 6 l/min with automatic compensation on forced ventilation modes and dedicated output
-----	--

---

Patient circuit	<ul style="list-style-type: none"> <li>• Double hose 150 cm. Adult/Paediatric patient circuit (expiratory valve on the ventilator)</li> <li>• Double hose 150 cm. Neonatal patient circuit (expiratory valve on the ventilator)</li> </ul>
-----------------	--

---

Software upgrade	USB 1 port
------------------	------------

## USER INTERFACE

Touch screen monitor	Module with TFT LED display with touch screen
----------------------	---

---

*Dimensions* 10,4"

*Displaying area* 262x163 mm

---

Display keyboard	<p>Keyboard for rapid access of functions. Encoder knob for:</p> <ul style="list-style-type: none"> <li>• selection, set up and confirmation of physiological breathing parameters</li> <li>• selection and direct activation of function</li> </ul>
------------------	--

---

Displaying and settings	<ul style="list-style-type: none"> <li>• Operative Mode setting</li> <li>• Visualization of alarm messages and signals</li> <li>• Setting and monitoring of physiological breathing parameters</li> <li>• Visualization of additional graphs and breathing parameters</li> <li>• MENU function for setting operation parameters</li> <li>• Activation of special functions</li> <li>• Visualization of operative mode, clock, date and time functions</li> <li>• Visualization of software version</li> </ul>
-------------------------	---

---

---

**Calibration Programs**

- Self Test
- Turbine Characterization
- Expiratory Flow Sensor Calibration
- Usage at High Altitude
- VTEc
- Nebulizer Enable
- ScreenShoot Enable
- Tourn Off

---

**MENU function - SETUP**

- Display (*Brightness, Energy Saving, Sound Volume, Touch Audio*)
- Date & Time
- Language
- Units (*Weight, Height, CO<sub>2</sub>, Pressure*)
- Default (*Erase Trend data, Erase Patient data, Setting & Ventilation Default*)
- Other (*NIV Enable, Power Failure, Apnea Time, Change Password, Save to USB*)
- IRMA/ISA (*Gas Sensor*)
- Supplementary Tests (*Expiratory Flow Sensor Calibration, O<sub>2</sub> Sensor Calibration*)
- Turn Off?
- Cancel

---

**MENU function - PATIENT DATA**

The PATIENT DATA can be set or deleted

---

**Alarm Limits**

PAW (cmH<sub>2</sub>O), PEEP (cmH<sub>2</sub>O), Vte (ml), VM (L/min), O<sub>2</sub> (%), RR (bpm), EtCO<sub>2</sub> (%)

---

**Displayed graphics**

- CURVES: Pressure (PAW) - Flow - Volume (Vte) - O<sub>2</sub> (CO<sub>2</sub> optional)
- LOOPS: Pressure / Volume - Flow / Volume - Pressure/Flow
- Graphics: INSP-EXP cycle
- Events
- Trends

---

*Events* Memory storage up to 100 machine events including the alarms.

---

*Trends* Storage capacity (72 h) of all measured parameters.

---

---

Physiological breathing parameters setting	Vti (ml), RR (bpm), I:E, Pause (%), PEEP (cmH <sub>2</sub> O), O <sub>2</sub> (%), Tr. I (L/min - cmH <sub>2</sub> O), SIGH (Sigh. Amp. (%), Sigh. Int. (b)), Vte (ml), PMax, Pmin, P <sub>insp</sub> (cmH <sub>2</sub> O), Slope, BACK-UP parameters, PS (cmH <sub>2</sub> O), RRsimv (bpm), Ti (s), Ti Max (s), Tr. E (%), CPAP (cmH <sub>2</sub> O), Pressure High - Low (cmH <sub>2</sub> O), Time High - Low (s).
--	--

---

- |                                     |  |
|-------------------------------------|--|
| <i>Range of measured parameters</i> | <ul style="list-style-type: none"> <li>• Respiratory rate (range: 0 ÷ 200 bpm)</li> <li>• I:E ratio (range 1:99 ÷ 99:1)</li> <li>• O<sub>2</sub>% (range: 0% ÷ 100%)</li> <li>• Tidal Volume: Vte, Vti (range: 0 ÷ 3000 ml)</li> <li>• Minute Volume (range: 0 ÷ 40 l/min)</li> <li>• PAW: peak, mean, plateau, PEEP (range -20 ÷ 80 cmH<sub>2</sub>O)</li> <li>• Inspiratory Peak Flow: Fi (range: 1 ÷ 190 l/min)</li> <li>• Expiratory Peak Flow: Fe (range: 1 ÷ 150 l/min)</li> <li>• T<sub>insp.</sub>, T<sub>exp</sub>, T<sub>pause</sub> (range 0.036 ÷ 10.9 sec)</li> <li>• Static and Dynamic compliance (range: 10 ÷ 150 ml/cmH<sub>2</sub>O)</li> <li>• Resistance (range: 0 ÷ 400 cmH<sub>2</sub>O/l/s)</li> <li>• EtCO<sub>2</sub>: with optional CO<sub>2</sub> module (range: 0 ÷ 10%)</li> <li>• Leak (%) (range: 0 ÷ 100%)</li> <li>• O<sub>2</sub> consumption (range: 0 ÷ 100l/min)</li> </ul> |
|-------------------------------------|--|
- 

<i>Displayed parameters</i>	PAW, PEEP, CPAP (cmH <sub>2</sub> O), RR (bpm), I:E, O <sub>2</sub> (% - l/min), Vte (ml), VM (L/min), EtCO <sub>2</sub> (%), MAP (cmH <sub>2</sub> O), P <sub>plateau</sub> (cmH <sub>2</sub> O), Fi, Fe (L/min), Ti, T <sub>pause</sub> , Te (sec.), Ri (cmH <sub>2</sub> O/l/s), Cs, Cd (ml/cmH <sub>2</sub> O), Leak ( % )
-----------------------------	--

---

<b>Flow sensor</b>	Magnetic perturbation (patented), reusable
--------------------	--

---

<i>Calibration</i>	Automatic (started by the operator)
--------------------	-------------------------------------

---

<i>Maintenance</i>	By steam or chemical disinfection
--------------------	-----------------------------------

---

<b>Oxymeter</b>	Electronic (value displayed in breathing parameters)
-----------------	--

---

<i>Calibration</i>	Automatic or started by the Operator
--------------------	--------------------------------------

---

<b>CO<sub>2</sub> analyzer</b>	Optional function (Sidestream or Mainstream module available)
--------------------------------	---

---

## ALARMS

Alarm types	<ul style="list-style-type: none"> <li>By MENU: with limits set by the operator</li> <li>By DEFAULT: the operator cannot set them up</li> </ul>
Alarm default setting	Present (Neonatal, Paediatric, Adult)
Alarm priority	High - Mean - Standby
Alarms visualization	Max 3 alarms simultaneously; additional alarms, scroll every 3-5 sec.

### *Alarms with limits set up by the operator*

Pressure of Airways	High – Low
Respiratory Rate	High – Low
Expiratory Volume	High – Low
Volume Minute	High – Low
PEEP	High – Low
O <sub>2</sub> Concentration	High – Low
EtCO <sub>2</sub>	High – Low (with optional CO <sub>2</sub> gas analyser)
On Battery	Alarm occurs in case of failure of external power supply
Apnoea	Low Rate (function of Apnoea BACK-UP)

### *System alarms*

Low Battery: 50% Remaining	Battery at 50%
Low Battery: 25% Remaining	Battery at 25%
Low Battery	10 Minutes
Battery Disconnected	Yes / No
Battery Overtemperature	Indication of exceeding the temperature limits inside the battery
Circuit Disconnected	Indication of patient circuit disconnected
O <sub>2</sub> Supply	Low (< 2,7 bar)
Turbine Failure	Signals in case of a blower fault condition
Turbine Overtemperature	Indication of exceeding the temperature limits inside the turbine
Turbine Overcurrent	Indication of exceeding the current limits inside the turbine
Maintenance	1000 hours
CO <sub>2</sub> Analyzer	Sampling Line Clogged, No Sampling Line, Replace Adapter, No Adapter, Unspecified Accuracy, Error, No Breaths, Low/High EtCO <sub>2</sub> .



## SELF-TEST alarms

Turbine	The correct functioning of the turbine is tested
Oxygen emptying	It is performed a washing of the remaining oxygen present within the lung ventilator, order to measure the offset of the oxygen sensor
INSP.- EXP. Flow sensor	Verification of EXP flow sensor operation
Pressure sensor	Verification of pressure sensor operation through control of PAW reading
Electrovalve	The correct functioning of electrovalve is tested
Patient circuit	Verification of patient circuit
Battery	Checking on battery power
Oxygen sensor	Cell condition
Acoustic alarm	Verification by the user of acoustic signal emission, the confirmation of the test is made by silencing of that alarm

## ACCESSORIES

### **Supplied Accessories**

- User's Manual
- Double hose patient circuit
- Antibacterial filter for patient circuit
- Nebulizer set
- Power cable
- Vehicular cable for 12 Vdc
- O<sub>2</sub> supply hose
- O<sub>2</sub> cell

### **Optional Accessories**

*Refer to price list*

SIARE applies the UNI EN ISO 13485:2016 Quality System and the 93/42/EEC.

**SIARE ENGINEERING INTERNATIONAL GROUP s.r.l.**

Via Pastore, 18 - Località Crespellano, 40053 Valsamoggia (BO), ITALY

Tel: +39 051 969802 - Fax: +39 051 969809

Email: [mail@siare.it](mailto:mail@siare.it) - <http://www.siare.it>